










1.0-1.2TM
PULL THRU
FEEL THE CLEAN


Single Patient Use Channel Cleaning Brush


1.0 -1.2 mm id

≈ 750-2200 mm

 Single use only
 USE WET
 Read Instructions before use
 NON STERILE




see www.pullthru.net.nz for more information


 F R Galantai Manufacturing Co. Ltd
 11 Bancroft Crescent, Glendene
 Auckland, 0602
 New Zealand
 ISO 13485:2016 Certified Site


 Advena Ltd.
 Tower Business Centre, 2nd Flr.,
 Tower Street, Swatar,
 BKR 4013 Malta

New Zealand Patent Number 540125,
 Other International Patents Pending
 Copyright © 2018 F R Galantai Manufacturing Co.Ltd.
 Trademarks are the property of their respective owners

